



Recombinant Human Zinc finger protein 610 (ZNF610)

Product Code	CSB-YP026890HU
Abbreviation	ZNF610
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9Z0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLCDEEAQKR KAKESGMALP QGRLTFMDVA IEF SQEEWKS LDPGQRALYR DVML ENYRNL VFLGICLPDL SIISMLKQRR EPLILQSQVK IVKNTDGREC VRSVNTGRSC VLG SNAENKP IKNQLGLTLE AHLSELQLFQ AGRKIYRSNQ VEKFTNHRSS VSPLQKISSS FTTHIFNKYR NDLIDFPLLP QEEKAYIRGK SYEYECSEDG EVFRVRASLT NHQVIHTAEK PYKCTECGKV FSRNSHLVEH WRIHTGQKPY KCSECDKVF N RNSNLARHQR IHTGEKPHKC NECGKAFREC SGLTTHLVIIH TGEKPYKCNE CGKNFRHKFS LTNHQRSHTA EKPYKCNECG KVFSLLSYLA RHQIIHSTEK PYKCNECGRA FHKRPGLMAH LLIHTGEKPY KCNECDKVFG RKLYLTNHQR IHTGERPYKC NACGKVF NQN PHLSRHRKIH AGENSLRTLQ ME
Source	Yeast
Target Names	ZNF610
Protein Names	Recommended name: Zinc finger protein 610
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.