



Recombinant Mouse WAS protein family homolog 1 (Wash1)

Product Code	CSB-BP025972MO
Abbreviation	Wash1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VDD8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MTPVKTKQCSL AGQLYAVPLI QPDLRREEAI QQVADALQYL QNISGDIFSR ISRRVELSRR QLQAISERSV LAQAKIEKIK GSKKAIKVFS SAKYPAPEHL QEYGSIFTGA LDPGLQRRPR YRIQSKHRPL DERALQEKLK YFPVCVNTKS EPEDEAEGL GGLPSNISSI SLLLLFNTTE NLYKKYVFLD PLAGAVTKTH TMLGTEEEKL FDAPLSISKR EQLERQAPEN YFYVPDLGQV PEIDVPSYLP DLPGVADDLM YSADLGGIA PSAPGAIPEL PAFHTEVAEP LQPELENEVL LAAPPPPPP PPPPPAPTA LVSTPQPPMF PDMATAAGQV AREEDSSSSM AHTASVQGAP KEVVDPPSSGR ATLLSIRQA GGIGKAKLRS VKERKLEKKK QKEQEVRAT SQGGDLMSDL FNKLVMRKKG ISGKGPSTGT SEGGGAFSR MSDSIPPLPP PQQPAGDEDE EDWES
Source	Baculovirus
Target Names	Washc1
Protein Names	Recommended name: WAS protein family homolog 1
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.