



# Recombinant Thymidine phosphorylase (deoA)

<b>Product Code</b>	CSB-EP025392EOD-B
<b>Abbreviation</b>	deoA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8XB35
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFLAQEIIRK KRDGHALSDE EIRFFINGIR DNTISEGQIA ALAMTIFFHD MTMPERVSLT MAMRDSGTVL DWKSLHLNGP IVDKHSTGGV GDVTSMLLGP MVAACGGYIP MISGRGLGHT GGTLDKLESI PGFDIFPDDN RFREIHKDVG VAIIGQTSSL APADKRFYAT RDITATVDSI PLITASILAK KLAEGLDALV MDVKVGSGAF MPTYELSEAL AEAIVGVANG AGVRTTALLT DMNQVLASSA GNAVEVREAV QFLTGEYRNP RLFDVTMALC VEMLISGKLA KDDAEARAKL QAVLDNGKAA EVFGRMVAAQ KGPTDFVENY AKYLPTAMLT KAVYADTEGF VSEMDTRALG MAVVAMGGGR RQASDTIDYS VGFTDMARLG DQVDGQRPLA VIHAKDENSW QEAAKAVKAA IKLADKAPES TPTVYRRISE
<b>Source</b>	E.coli
<b>Target Names</b>	deoA
<b>Protein Names</b>	Recommended name: Thymidine phosphorylase EC= 2.4.2.4 Alternative name(s): TdRPase
<b>Expression Region</b>	1-440
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.