



Recombinant Mouse Ras and Rab interactor 1 (Rin1)

| | |
|--------------------------|--|
| Product Code | CSB-BP019726MO |
| Abbreviation | Rin1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q921Q7 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MEDPGETGAH PLGATNLNFV PGHQQKEKPS TDPLYDTPDT RGVQAGGSQQ PARTVSLRER LLITRPVWLQ LRANAAAALH VLRTEPPGTF LVRKSNTRQC QALCVRLPEA SGPSFVSSHY IEESTGGVSL EGSELMFQDL VQLICGYCRT RAIHQAATHK ELEAISHLGM EFWSSSLNTK DQQRPSEAPP IPRLKARSPQ ELDQGTGAAL CFFNPLFPGD LGPTKREKFK RSFKVRVSTE TSSPLSPPAV PPPVPVLPG TSSSQTERLP PRQLLQRESS VGYRVPGSAA SPCLPPLPSL QEVDCSPSS SEEEGSSGSP TTSPRLSRPR HRRPLLRSMS SAFCSLLAPE RQVGRAATML MQNRYTAVGQ LVQDLLTQVR AGPEPRELQG IRQALSRARA MLSAELGPEK LLPPERLELV LEKSLHRSVL KPLRPILAAR LRRRLSADGS LGRLAEGFRL ARTQGGPAGF SHLTLSPPVE TEQVRQKLLQ LLRAYSPSAQ VKWLLQACKL LYTALKSQAG ENAGADEFLP LLSLVLAQCD LPDLLLEAEY MSELLEPTLL TGEAGGYLTS LSASLALLSG LSQARALPLS PAQELQRSLA LWEQRRLPAT HSFQHLLRVA YQDPSTGCTS KTLAVPPGSS IATLSQLCAT KFRVTQPDAF GLFLYKDQGY HRLPPEALAH RLPATGYLIY RRAERPETQG AVAEKAKTGS KGPEAGAWEE ETGGLNREGK PRIAVDQEGK DQARGGHIGP EEQKAEGSQA LEE |
| Source | Baculovirus |
| Target Names | Rin1 |
| Protein Names | Recommended name: Ras and Rab interactor 1 Alternative name(s): Ras interaction/interference protein 1 |
| Expression Region | 1-763 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.