



Recombinant Human Zinc finger protein 614 (ZNF614)

Product Code	CSB-MP026893HU
Abbreviation	ZNF614
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N883
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MIKTQESLTL EDVAVEFSWE EWQLLDTAQK NLYRDVMVEN YNHLVSLGYQ TSKPDVLSKL AHGQEPWTTD AKIQNKNC PG IGKVDSHLQE HSPNQRL LKS VQQCNGQNTL RNIVHLSKTH FPIVQNHDTF DLYRKNL KSS LSLINQKRRH GINNPVEFIG GEKTL LHGKH ERTHTKTRFS ENAKCIHTKF QVFKHQ RTQK IEKPHACIEC EQTFLRKS QL IYHENICIQE NPGSGQCEKL SRSVLFTKHL KTNTTDKICI PNEYRKGSTV KSSLITHQQT HTEEKSYMCS ECGKGFTMKR YLIAHQ RTHS GEKPYVCKEC GKGFTVKS NL IVHQ RTHTGE KPYICSECGK GFTMKRYLVV HQRTHTGEKP YMCSECGKGF TVKSNLIVHQ RSHTGEKSYI CSECGKGFTV KRTLVIHQ RT HTGEKSYICN ECGKGFTTKR TLIHQ RTHT GEKPYECNEC GKAFSQKICL IQHERCHTGK TPFVCTECGK SYSHKYGLIT HQRIHTGEKP YECNECGKAF TTKSVLNVHQ RTHTGERPYG CSDCEKAFSH LSNLVKHKKM HTREMGRISQ VENSNGESQ LLPYK
Source	Mammalian cell
Target Names	ZNF614
Protein Names	Recommended name: Zinc finger protein 614
Expression Region	1-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.