



# Recombinant Human Fibroblast growth factor 14 (FGF14)

<b>Product Code</b>	CSB-MP008620HU
<b>Abbreviation</b>	FGF14
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q92915
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVKPVPLFRRTDFKLLLCNHKDLFFLRVSKLLDCFSPKSMWFLWNIFSKGTHM LQCLCGKSLKKNKNPTDPQLKGIVTRLYCRQGYLQMHDPDGDGTDKDDSTN STLFNLIPVGLRVVAIQGVKTGLYIAMNGEGYLYPSELFTEPECKFKESVFENYYV IYSSMLYRQQESGRAWFLGLNKEGQAMKGNRVKKTTPAAHFLPKPLEVAMYR EPSLHDVGETVPKPGVTPSKSTSASAIMNGGKPVNKSSTT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FGF14
<b>Protein Names</b>	Recommended name: Fibroblast growth factor 14 Short name= FGF-14Alternative name(s): Fibroblast growth factor homologous factor 4 Short name= FHF-4
<b>Expression Region</b>	1-252
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length of Isoform 2
<b>Target Details</b>	This protein is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. A mutation in this gene is associated with autosomal dominant cerebral ataxia. Alternatively spliced transcript variants have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.