



Recombinant Human Zinc finger protein 830 (ZNF830)

Product Code	CSB-MP027073HU
Abbreviation	ZNF830
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NB3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ASSASARTP AGKRVINQEE LRRLMKEKQR LSTSRKRIES PFAKYNRLGQ LSCALCNTPV KSELLWQTHV LGKQHREKVA ELKGAKEASQ GSSASSAPHS VKRKAPDADD QDVKRAKATL VPQVQPSTSA WTTNFDKIGK EFIRATPSKP SGLSLLPDYE DEEEEEEEEE GDGERKRGDA SKPLSDAQGK EHSVSSSREV TSSVLPNDFF STNPPKAPII PHSGSIEKAE IHEKVVERRE NTAEALPEGF FDDPEVDARV RKVDAPKDQM DKEWDEFQKA MRQVNTISEA IVAEEDEEGR LDRQIGEIDE QIECYRRVEK LRNRQDEIKN KLKEILTIKE LQKKEENAD SDDEGELQDL LSQDWRVKGA LL
Source	Mammalian cell
Target Names	ZNF830
Protein Names	Recommended name: Zinc finger protein 830 Alternative name(s): Coiled-coil domain-containing protein 16
Expression Region	2-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.