



# Recombinant Human Mitogen-activated protein kinase kinase kinase 3 (MAP3K3)

<b>Product Code</b>	CSB-EP013425HU
<b>Abbreviation</b>	MAP3K3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q99759
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDEQEALNSI MNDLVALQMN RRHRMPGYET MKNKDTGHSN RQSDVRIKFE HNGERRIIAF SRPVKYEDVE HKVTTVFGQP LDLHYMNNEL SILLKNQDDL DKAIDILDRS SSMKSLRILL LSQDRNHNS SPHSGVSRQV RIKASQSAGD INTIYQPPEP RSRHLSVSSQ NPGRSSPPPG YVPERQQHIA RQGSYTSINS EGEFIPETSE QCMLDPLSSA ENSLSGSCQS LDRSADSPSF RKSRRMSRAQS FPDNRQEYSD RETQLYDKGV KGGTYPRRYH VSVHHKDYSYD GRRTFPRIRR HQQNLFTLVP SSRSLSTNGE NMGLAVQYLD PRGRLRSADS ENALSVQERN VPTKSPSAPI NWRRGKLLGQ GAFGRVYLCY DVDTGRELAS KQVQFDPDSP ETSKEVSALE CEIQLLKNLQ HERIVQYYGC LRDRAEKTTLT IFMEYMPGGS VKDQLKAYGA LTESVTRKYT RQILEGMSYL HSNMIVHRDI KGANILRDSA GNVKLGDFGA SKRLQTICMS GTGMRSVTGT PYWMSPEVIS GEGYGRKADV WSLGCTVVEM LTEKPPWAEY EAMAAIFKIA TQPTNPQLPS HISEHGRDFL RRIFVEARQR PSAEELLTHH FAQLMY
<b>Source</b>	E.coli
<b>Target Names</b>	MAP3K3
<b>Protein Names</b>	Recommended name: Mitogen-activated protein kinase kinase kinase 3 EC=2.7.11.25 Alternative name(s): MAPK/ERK kinase kinase 3 Short name= MEK kinase 3 Short name= MEKK 3
<b>Expression Region</b>	1-626
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene product is a 626-amino acid polypeptide that is 96.5% identical to mouse Mekk3. Its catalytic domain is closely related to those of several other kinases, including mouse Mekk2, tobacco NPK, and yeast Ste11. Northern blot analysis revealed a 4.6-kb transcript that appears to be ubiquitously expressed.



This protein directly regulates the stress-activated protein kinase (SAPK) and extracellular signal-regulated protein kinase (ERK) pathways by activating SEK and MEK1/2 respectively; it does not regulate the p38 pathway. In cotransfection assays, it enhanced transcription from a nuclear factor kappa-B (NFkB)-dependent reporter gene, consistent with a role in the SAPK pathway. Alternatively spliced transcript variants encoding distinct isoforms have been observed.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.