



Recombinant Human Zinc finger protein 238 (ZNF238)

Product Code	CSB-MP026604HU
Abbreviation	ZNF238
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99592
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEFPDHSRHL LQCLSEQRHQ GFLCDCTVLV GDAQFRAHRA VLASCSMYFH LFYKDQLDKR DIVHLNSDIV TAPAFALLLE FMYEGKLQFK DLPIEDVLAA ASYLHMYDIV KVCKKKLKEK ATTEADSTKK EEDASSCSKD VESLSDGSSH IAGDLPSDED EGEDEKLNIL PSKRDLAAEP GNMWMRLPSD SAGIPQAGGE AEPHATAAGK TVASPCSSTE SLSQRSVTSV RDSADVDCVL DLSVKSSLSG VENLNSSYFS SQDVLRSNLV QVKVEKEASC DESDVGNDY DMEHSTVKES VSTNNRVQYE PAHLAPLRED SVLRELDRED KASDDEMTP ESERVQVEGG MESSLLPYVS NILSPAGQIF MCPLCNKVFP SPHILQIHLS THFREQDGIR SKPAADVNP TCSLCGKTFS CMYTLKRHER THSGEKPYTC TQCGKSFQYS HNLRSRHAVVH TREKPHACKW CERRFTQSGD LYRHIRKFHC ELVNSLSVKS EALSPTVRD WTLEDSSQEL WK
Source	Mammalian cell
Target Names	ZBTB18
Protein Names	Recommended name: Zinc finger protein 238 Alternative name(s): 58 kDa repressor protein Transcriptional repressor RP58 Translin-associated zinc finger protein 1 Short name= TAZ-1 Zinc finger and BTB domain-containing protein 18
Expression Region	1-522
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.