



Recombinant Callithrix jacchus Deleted in azoospermia-like (DAZL)

Product Code	CSB-MP006514CYL
Abbreviation	DAZL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BGN8
Product Type	Recombinant Protein
Immunogen Species	Callithrix jacchus (White-tufted-ear marmoset)
Purity	>85% (SDS-PAGE)
Sequence	MSAANPETPN STISREANTQ SSSAATSQGY VLPEGKIMPN TVFVGGIDVR MDETEIRGFF ARYGSVKEVK IITDRTGVSK GYGFVSFFND VDVQKIVESQ INFHGKCLKL GPAIRKQNLK AYHVQPRPLV FNHPPPPQFQ NVWSNPNTET YMHPPTTMNP VTQYVQAYPP YPNSPVQVIT GYQLPVYNYQ MPPQWPVGEQ RSYVVPPAYS SVNYCNEID PGAEVVNEC SVCEATPPSG NGPQKKSVDK SIQTVVSLF NPENSLRNS VVTQDDYFRD KRVHHFRRSR AVLKSV
Source	Mammalian cell
Target Names	DAZL
Protein Names	Recommended name: Deleted in azoospermia-like Alternative name(s): DAZ-like autosomal Deleted in azoospermia-like 1
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.