



# Recombinant Human Histone H4 transcription factor (HINFP)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-BP010361HU  |
| <b>Abbreviation</b>      | HINFP   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9BQA5  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MPPPGKVPRK ENLWLQCEWG SCSFVCSTME KFFEHVVTQHL<br>QQHLHGSGEE EEEEEEDDPL EEEFSCWQEQ CGFCSLDSSA DLIRHVFYFHC<br>YHTKLKQWGL QALQSQADLG PCILDFQSRN VIPDIPDHFL CLWEHCENSF<br>DNPEWFYRHV EAHSLCCEYE AVGKDNPVVL CGWKGCTCTF KDRSKLREHL<br>RSHTQEKVVA CPTCGGMFAN NTKFLDHIRR QTSLDQQHFQ CSHCSKRFAT<br>ERLLRDHMRN HVNHYKCPLC DMTCLPSSL RNHMRFRHSE<br>DRPFKDCDD YSCKNLIDLQ KHLDTHTSEEP AYRCDFENCT FSARSLCSIK<br>SHYRKVHEGD SEPRYKCHVC DKCFTRGNL TVHLRKKHQF<br>KWPSGHPRFR YKEHEDGYMR LQLVRYESVE LTQQLLRQPQ EGSGLGTSLN<br>ESSLQGIILE TVPGEPGRKE EEEEGKQSEG TALSASQDNP SSVIHVVNQT<br>NAQQQQEIVY YVLSEAPGEP PPAPEPPSGG IMEKLQGGIAE EPEIQMV |
| <b>Source</b>            | Baculovirus   |
| <b>Target Names</b>      | HINFP   |
| <b>Protein Names</b>     | Recommended name: Histone H4 transcription factor Alternative name(s):<br>Histone nuclear factor P Short name= HiNF-P MBD2-interacting zinc finger<br>protein Methyl-CpG-binding protein 2-interacting zinc finger protein  |
| <b>Expression Region</b> | 1-517   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.