



Recombinant Human Histone H4 transcription factor (HINFP)

Product Code	CSB-BP010361HU
Abbreviation	HINFP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BQA5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPPPGKVPRK ENLWLQCEWG SCSFVCSTME KFFEHVVTQHL QQHLHGSGEE EEEEEEDDPL EEEFSLWQE CGFCSLDSSA DLIRHVYFHC YHTKLKQWGL QALQSQADLG PCILDFQSRN VIPDIPDHFL CLWEHCENSF DNPEWFYRHV EAHSLCCEYE AVGKDNPVVL CGWKGCTCTF KDRSKLREHL RSHTQEKVVA CPTCGGMFAN NTKFLDHIRR QTSLDQQHFQ CSHCSKRFAT ERLLRDHMRN HVNHYKCPLC DMTCLPSSL RNHMRFRHSE DRPFKCDCCD YSCKNLIDLQ KHLDTHEEP AYRCDFENCT FSARSLCSIK SHYRKVHEGD SEPRYKCHVC DKCFTRGNL TVHLRKKHQF KWPSGHPRFR YKEHEDGYMR LQLVRYESVE LTQQLLRQPQ EGSGLGTSLN ESSLQGIILE TVPGEPGRKE EEEEGKGSEG TALSASQDNP SSVIHVVNQT NAQQQQEIVY YVLSEAPGEP PPAPEPPSGG IMEKLQGIAE EPEIQMV
Source	Baculovirus
Target Names	HINFP
Protein Names	Recommended name: Histone H4 transcription factor Alternative name(s): Histone nuclear factor P Short name= HiNF-P MBD2-interacting zinc finger protein Methyl-CpG-binding protein 2-interacting zinc finger protein
Expression Region	1-517
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.