



# Recombinant Bovine Allograft inflammatory factor 1 (AIF1)

<b>Product Code</b>	CSB-MP001490BO
<b>Abbreviation</b>	AIF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BDK2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SETRDLQGG KAFGLRKAQQ EERINEINQQ FLDDPKYSSD EDLPSKLEAF KKKYMFDLN EDGGIDIMSL KRMMEKLGVP KTHLELKLI MEVSSGPGET FSYSDFLKMM LGKRSAILKM ILMYEEKARE QEKPTGLPAK KAISELP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	AIF1
<b>Protein Names</b>	Recommended name: Allograft inflammatory factor 1 Short name= AIF-1
<b>Expression Region</b>	2-147
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene is induced by cytokines and interferon. Its protein product is thought to be involved in negative regulation of growth of vascular smooth muscle cells, which contributes to the anti-inflammatory response to vessel wall trauma. Three transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.