



Recombinant Human Activating signal cointegrator 1 complex subunit 2 (ASCC2)

Product Code	CSB-MP002197HU
Abbreviation	ASCC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H118
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPALPLDQLQ ITHKDPKTGK LRTSPALHPE QKADRYFVLY KPPPKNIPALVEEYLERAT FVANDLDWLL ALPHDKFWCQ VIFDETLQKC LDSYLRYVPRKFDEGVASAP EVVDMQKRLH RSVFLTFLRM STHKESKDHFISSAFGEILYNNFLFDIPK ILDLCVLF GK GNSPLLQKMI GNIFTQQPSY YSDLDELPTILQVFSNILQ HCGLQGDGAN TTPQKLEERG RLTPSDMPLL ELKDIVLYLCDTCTTLWAFL DIFPLACQTF QKHDFCYRLA SFYEAAIPEM ESAIKKRRLEDSKLLGDLWQ RLSHSRKKLM EIFHIILNQI CLLPILESSC DNIQGFIEEF LQIFSSLLQE KRFLRDYDAL FPVAEDISLL QQASSVLDLRTAYILQAVE SAWEGVDRRK ATDAKDPSVI EEPNGEPNGV TVTAEAVSQA SSHPENSEEE ECMGAAAAGV PAMCGVELDS LISQVKDLLP DLGEGFILAC LEYYHYDPEQVINNILEERL APTLSQLDRN LDREMCPDPT PLLTSRHNVF QNDEFDVFSDSVDLSRVHK GKSTRKEENT RSLNDKRAV AAQRQRYEQY SVVVEEVPLQPGESLPYHSV YYEDEYDDTY DGNQVGANDA DSDDELISRR PFTIPQVLRTKVPREGQEED DDDEEDDADE EAPKPDHFVQ DPAVLREKAE ARRMAFLAKKGYRHDSSTAV AGSPRGHGQS RETTQERRKK EANKATRANH NRRTMADRKR SKGMIPS
Source	Mammalian cell
Target Names	ASCC2
Protein Names	Recommended name: Activating signal cointegrator 1 complex subunit 2 Alternative name(s): ASC-1 complex subunit p100 Trip4 complex subunit p100
Expression Region	1-757
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.