



# Recombinant Human Activating signal cointegrator 1 complex subunit 2 (ASCC2)

<b>Product Code</b>	CSB-MP002197HU
<b>Abbreviation</b>	ASCC2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H118
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPALPLDQLQ ITHKDPKTGK LRTSPALHPE QKADRYFVLY KPPPKNIPALVEEYLERAT FVANDLDWLL ALPHDKFWCQ VIFDETLQKC LDSYLRYVPRKFDEGVASAP EVVDMQKRLH RSVFLTFLRM STHKESKDFH ISPSAFGEILYNNFLFDIPK ILDLCVLF GK GNSPLLQKMI GNIFTQQPSY YSDLDETLPTILQVFSNILQ HCGLQGDGAN TTPQKLEERG RLTPSDMPLL ELKDIVLYLCDTCTTLWAFL DIFPLACQTF QKHDFCYRLA SFYEAAIPEM ESAIKKRRLED SKLLGDLWQ RLSHSRKKLM EIFHIILNQI CLLPILESSC DNIQGFIEEF LQIFSSLLQE KRFLRDYDAL FPVAEDISLL QQASSVLDET RTAYILQAVE SAWEGVDRRK ATDAKDPSVI EEPNGEPNGV TVTAEAVSQA SSHPENSEEE ECMGAAAAGV PAMCGVELDS LISQVKDLLP DLGEGFILAC LEYYHYDPEQ VINNILEERL APTLSQLDRN LDREMCPDPT PLLTSRHNVF QNDEFDVFSDSVDLRVHK GKSTRKEENT RSLNDKRAV AAQRQRYEQY SVVVEEVPLQPGESLPYHSV YYEDEYDDTY DGNQVGANDA DSDELISRR PFTIPQVLRTKVPREGQEED DDDEEDDADE EAPKPDHFVQ DPAVLREKAE ARRMAFLAKK GYRHDSSTAV AGSPRGHGQS RETTQERRKK EANKATRANH NRRTMADRKR SKGMIPS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ASCC2
<b>Protein Names</b>	Recommended name: Activating signal cointegrator 1 complex subunit 2 Alternative name(s): ASC-1 complex subunit p100 Trip4 complex subunit p100
<b>Expression Region</b>	1-757
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.