



# Recombinant Arabidopsis thaliana Mitogen-activated protein kinase 10 (MPK10)

<b>Product Code</b>	CSB-EP013449DOA
<b>Abbreviation</b>	MPK10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9M1Z5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEPTNDAETL ETQGEVTTAI WPSSQILKTT IDIPGTLSHD GRYIQYNLFG HIFELPAKYK PPIRPIGRGA CGIVCSAVDS ETNEKVAIKK ITQVFDNTIE AKRTLREIKL LRHFDHENIV AIRDVILPPQ RDSFEDVYIV NELMEFDLYR TLKSDQELTK DHGMYFMYQI LRGLKYIHSA NVLHRDLKPS NLLLSTQCDL KICDFGLARA TPESNLMTEY VVTRWYRAPE LLLGSSDYTA AIDVWSVGC FMEIMNREPL FPGKDQVNQL RLLLELIGTP SEEELGSLSE YAKRYIRQLP TLPRQSFTEK FPNVPLAID LVEKMLTFDP KQRISVKEAL AHPYLSSFHD ITDEPECSEP FNFDLDEHPF SEEQFRELIY CEALAFNPET SND
<b>Source</b>	E.coli
<b>Target Names</b>	MPK10
<b>Protein Names</b>	Recommended name: Mitogen-activated protein kinase 10 Short name= AtMPK10 Short name= MAP kinase 10 EC= 2.7.11.24
<b>Expression Region</b>	1-393
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.