



Recombinant Human Zinc finger protein 2 (ZNF2)

Product Code	CSB-EP026572HU-B
Abbreviation	ZNF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSG1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAVSPTTRC QESVTFEDVA VVFTDEEWSR LVPIQRDLYK EVMLENYNSI VSLGLPVPQP DVIFQLKRGD KPWMVDLHGS EEREWPESVS LDWETKPEIH DASDKKSEGS LRECLGRQSP LCPKFEVHTP NGRMGTEKQS PSGETRKKSL SRDKGLRRRR SALSREILTK ERHQECSDCG KTFFDHSSLT RHQRHTTGEK PYDCRECGKA FSHRSSLSRH LMSHTGESPY ECSVCSKAFF DRSSLTVHQR IHTGEKPFQC NECGKAFFDR SSLTRHQRIH TGESPYECHQ CGKAFSQKSI LTRHQLIHTG RKPYECNECG KAFYGVSSLN RHQKAHAGDP RYQCNECGKA FFDRSSLTQH QKIHTGDKPY ECSECGKAFFS QRCRLTRHQR VHTGEKPFEC TVCGKVFSSK SSVIQHQRRY AKQGID
Source	E.coli
Target Names	ZNF2
Protein Names	Recommended name: Zinc finger protein 2 Alternative name(s): Zinc finger protein 2.2 Zinc finger protein 661
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.