



Recombinant Chicken Secreted frizzled-related protein 1 (SFRP1)

Product Code	CSB-EP021138CH-B
Abbreviation	SFRP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9DEQ4
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	SEYDYVSYQ SDLGPYPGGR FYTKPHQCVA IPADLRLCHS VGYDKMVLPN LLDHETMAEV KHQASSWVPL LNKNCHEMGTO VFLCSLFAPV CLDRPVYPCR WLCEAVRDSC EPVMQFFGFF WPEMLKCDQF PQDYVCIAMT TPNATEVSRP KGTTCPPCD NEMKSEAIVE HLCASEFALK MTKIEVKKEN GDKVIIPRKR KALKLGPVK KNLKLVLLL KNGADCPCHQ LDNLGHHFLI MGRQVKTQHL LTAIYKWDKK NKEFKKFMKK VKAPDCPTFP SVFK
Source	E.coli
Target Names	SFRP1
Protein Names	Recommended name: Secreted frizzled-related protein 1 Short name= CsFRP1 Short name= sFRP-1
Expression Region	32-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. Members of this family act as soluble modulators of Wnt signaling; epigenetic silencing of SFRP genes leads to deregulated activation of the Wnt-pathway which is associated with cancer. This gene may also be involved in determining the polarity of photoreceptor cells in the retina.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.