



Recombinant Human RING finger protein 39 (RNF39)

Product Code	CSB-BP019888HU
Abbreviation	RNF39
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H2S5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MWWRDLTRLR LWLKREAIPG EGRKAAKVNA GVGEKGIYTA SSRGGPPSAR SKAVTVVAEG AASRSWLSMD APELGPGLVE RLEQLATCPL CGGSFEDPVL LACEHSFCRA CLARRWGTPP ATGTEASPTA CPCCGLPCPR RSLRSNVRLA VEVRI SREL R EKLAEPGARA GRRRGRIPT MGCLDLPGED MRKTWRRFEV PTSKSSNSED DLPEDYPVVK KMLHRLTADL TLDPGTAHRR LLISADRRSV QLAPPGTPAP PDGPKRFDQL PAVLGAQQFG AGRHCWEVET ADAASCRDSS GEDADDEESH YAVGAAGESV QRKGCVRLCP AGAVWAVEGR GGRLWALTAP EPTLLGGVEP PPRRIRVDLD WERGRVAFYD GRSLDLLYAF QAPGPLGERI FPLFCTCDPR APLRIVPAES
Source	Baculovirus
Target Names	RNF39
Protein Names	Recommended name: RING finger protein 39 Alternative name(s): Protein HZFW
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene lies within the major histocompatibility complex class I region on chromosome 6. Studies of a similar rat protein suggest that this gene encodes a protein that plays a role in an early phase of synaptic plasticity. Multiple transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.