



# Recombinant Human RING finger protein 39 (RNF39)

<b>Product Code</b>	CSB-MP019888HU
<b>Abbreviation</b>	RNF39
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H2S5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MWWRDLTRLR LWLKREAIPG EGRKAAKVNA GVGEKGIYTA SSRGGPPSAR SKAVTVVAEG AASRSWLSMD APELGPLVE RLEQLATCPL CGGSFEDPVL LACEHSFCRA CLARRWGTPP ATGTEASPTA CPCCGLPCPR RSLRSNVRLA VEVRI SREL R EKLAEPGARA GRRRGGRIPT MGCLDLPGED MRKTWRRFEV PTSKSSNSE D LPEDYPVVK KMLHRLTADL TLDPGTAHRR LLISADRRSV QLAPPGTPAP PDGPKRFDQL PAVLGAQGFAGRHCWEVET ADAASCRDSS GEDADDEESH YAVGAAGESV QRKGCVRLCP AGAVWAVEGR GGRLWALTAP EPTLLGGVEP PRRIRVDLD WERGRVAFYD GRSLDLLYAF QAPGPLGERI FPLFCTCDPR APLRIVPAES
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RNF39
<b>Protein Names</b>	Recommended name: RING finger protein 39 Alternative name(s): Protein HZFW
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene lies within the major histocompatibility complex class I region on chromosome 6. Studies of a similar rat protein suggest that this gene encodes a protein that plays a role in an early phase of synaptic plasticity. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.