



Recombinant Human Zinc finger protein 408 (ZNF408)

Product Code	CSB-MP026728HU
Abbreviation	ZNF408
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H9D4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEAEELLLE GKKALQLARE PRLGLDLGWN PSGEGCTQGL KDVPPEPTRD ILALKSLPRG LALGPSLAKE QRLGVWCVDG PLQPGLLWGP LEEESASKEK GEGVKPRQEE NLSLGPWGDV CACEQSSGWT SLVQRGRLES EGNVAPVRIS ERLHLQVYQL VLPGSELLLW PQPSSEGPSL TQPGLDKEAA VAVVTEVESA VQQEVASPGE DAAEPCIDPG SQSPSGIQAE NMVSPGLKFP TQDRISKDSQ PLGPLLQDGD VDEECPAQAQ MPPELQSNSA TQQDPDGSGA SFSSSARGTQ PHGYLAKKLH SPSDQCPPRA KTPEPGAQQS GFPTLSRSPP GPAGSSPKQG RRYRCGECGK AFLQLCHLKK HAFVHTGHKP FLCTECGKSY SSEESFKAHM LGHRGVRPFP CPQCDKAYGT QRDLKEHQVV HSGARPFACD QCGKAFARRP SLRLHRKTHQ VPAAPAPCPC PVCGRPLANQ GSLRNHMRLH TGEKPFLLCPH CGRAFRQRGN LRGHLRLHTG ERPYRCPHCA DAFPQLPELR RHLISHTGEA HLCPCVCGKAL RDPHTLRAHE RLHSGERPFP CPQCGRAYTL ATKLRRHLKS HLEDKPYRCP TCGMGYTLPQ SLRRHQLSHR PEAPCSPPSV PSAASEPTVV LLQAEPQLLD THREEEVSPA RDVVEVTISE SKEKCFVVPE EPDAAPSLVL IHKDMGLGAW AEVVEVEMGT
Source	Mammalian cell
Target Names	ZNF408
Protein Names	Recommended name: Zinc finger protein 408 Alternative name(s): PR domain zinc finger protein 17
Expression Region	1-720
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.