



# Recombinant Human Zinc finger protein 408 (ZNF408)

<b>Product Code</b>	CSB-YP026728HU
<b>Abbreviation</b>	ZNF408
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H9D4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEEAEELLLE GKKALQLARE PRLGLDLGWN PSGEGCTQGL KDVPPEPTRD ILALKSLPRG LALGPSLAKE QRLGVWCVDG PLQPGLLWGP LEEESASKEK GEGVKPRQEE NLSLGPWGDV CACEQSSGWT SLVQRGRLES EGNVAPVRIS ERLHLQVYQL VLPGSELLLW PQPSSEGPSL TQPGLDKEAA VAVVTEVESA VQQEVASPGE DAAEPCIDPG SQSPSGIQAE NMVSPGLKFP TQDRISKDSQ PLGPLLQDGD VDEECPAQAQ MPPELQSNSA TQQDPDGSGA SFSSSARGTQ PHGYLAKKLH SPSDQCPPRA KTPEPGAQQS GFPTLSRSPP GPAGSSPKQG RRYRCGECGK AFLQLCHLKK HAFVHTGHKP FLCTECGKSY SSEESFKAHM LGHRGVRPFP CPQCDKAYGT QRDLKEHQVV HSGARPFACD QCGKAFARRP SLRLHRKTHQ VPAAPAPCPC PVCGRPLANQ GSLRNHMRLH TGEKPFLLCPH CGRAFRQRGN LRGHLRLHTG ERPYRCPHCA DAFPQLPELR RHLISHTGEA HLCPCVCGKAL RDPHTLRAHE RLHSGERPFP CPQCGRAYTL ATKLRRHLKS HLEDKPYRCP TCGMGYTLPQ SLRRHQLSHR PEAPCSPPSV PSAASEPTVV LLQAEPQLLD THREEEVSPA RDVVEVTISE SKEKCFVVPE EPDAAPSLVL IHKDMGLGAW AEVVEVEMGT
<b>Source</b>	Yeast
<b>Target Names</b>	ZNF408
<b>Protein Names</b>	Recommended name: Zinc finger protein 408 Alternative name(s): PR domain zinc finger protein 17
<b>Expression Region</b>	1-720
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.