



# Recombinant Human Myoneurin (MYNN)

<b>Product Code</b>	CSB-EP015330HU
<b>Abbreviation</b>	MYNN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NPC7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MQYSHHCEHL LERLNKQREA GFLCDCTIVI GEFQFKAHRN VLASFSEYFG            AIYRSTSENN VFLDQSQVKA DGFQKLLEFI YTGTLNLDSW NVKEIHQAAD            YLKVEEVVTK CKIKMEDFAF IANPSSTEIS SITGNIELNQ QTCLLTLRDY            NNREKSEVST DLIQANPKQG ALAKKSSQTK KKKKAFNSPK TGQNKTQVQYP            SDILENASVE LFLDANKLPT PVVEQVAQIN DNSELELTSV VENTFPAQDI            VHTVTVKRKR GKSQPNCALK EHSMSNIASV KSPYEAENSG EELDQRYSKA            KPMCNTCGKV FSEASSLRRH MRIHKGVPY VCHLCGKAFT QCNQLKTHVR            THTGEKPYKC ELCDKGFAQK CQLVFHSRMH HGEEKPYKCD VCNLQFATSS            NLKIHARKHS GEKPYVCDRC GQRFAQASTL TYHVRRTGE KPYVCDTCGK            AFAVSSSLIT HSRKHTGEKP YICGICGKSF ISSGELNKH FRSHTGERPFI            CELCGNSYTD IKNLKHKTK VHSGADKTLD SSAEDHTLSE QDSIQKSPLS            ETMDVKPSDM TLPLALPLGT EDHHMLLPVT DTQSPTSDDL LRSTVNGYSE            PQLIFLQQLY         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	MYNN
<b>Protein Names</b>	Recommended name: Myoneurin Alternative name(s): Zinc finger and BTB domain-containing protein 31
<b>Expression Region</b>	1-610
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.