



Recombinant Human Myoneurin (MYNN)

Product Code	CSB-MP015330HU
Abbreviation	MYNN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NPC7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQYSHHCEHL LERLNKQREA GFLCDCTIVI GEFQFKAHRN VLASFSEYFG AIYRSTSENN VFLDQSQVKA DGFQKLLEFI YTGTLNLDSW NVKEIHQAAD YLKVEEVVTK CKIKMEDFAF IANPSSTEIS SITGNIELNQ QTCLLTRDY NNREKSEVST DLIQANPKQG ALAKKSSQTK KKKKAFNSPK TGQNKTQVQYP SDILENASVE LFLDANKLPT PVVEQVAQIN DNSELELTSV VENTFPAQDI VHTVTVKRKR GKSQPNCALK EHSMSNIASV KSPYEAENSG EELDQRYSKA KPMCNTCGKV FSEASSLRRH MRIHKGVPY VCHLCGKAFT QCNQLKTHVR THTGEKPYKC ELCDKGFAQK CQLVFHSRMH HGEEKPYKCD VCNLQFATSS NLKIHARKHS GEKPYVCDRC GQRFAQASTL TYHVRRTGE KPYVCDTCGK AFAVSSSLIT HSRKHTGEKP YICGICGKSF ISSGELNKH FRSHTGERPFI CELCGNSYTD IKNLKKHKTK VHSGADKTLD SSAEDHTLSE QDSIQKSPLS ETMDVKPSDM TLPLALPLGT EDHHMLLPVT DTQSPTS DTL LRSTVNGYSE PQLIFLQQLY</p>
Source	Mammalian cell
Target Names	MYNN
Protein Names	Recommended name: Myoneurin Alternative name(s): Zinc finger and BTB domain-containing protein 31
Expression Region	1-610
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.