



Recombinant Human Cysteine-rich PDZ-binding protein (CRIPT)

Product Code	CSB-MP005973HU
Abbreviation	CRIPT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P021
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVCEKCEKKL GTVITPDTWK DGARNTTESG GRKLNENKAL TSKKARFDPY GKNKFSTCRI CKSSVHQP GS HYCQG CAYKK GICAMCGKKV LDTKNYKQTS V
Source	Mammalian cell
Target Names	CRIPT
Protein Names	Recommended name: Cysteine-rich PDZ-binding protein Alternative name(s): Cysteine-rich interactor of PDZ three Short name= Cysteine-rich interactor of PDZ3
Expression Region	1-101
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.