



Recombinant Rat Inositol hexakisphosphate kinase 2 (Ip6k2)

Product Code	CSB-EP011771RA-B
Abbreviation	Ip6k2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9R0U1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MSPAFRTMDV EPRTKGILLE PFVHQVGGHS CVLRFNETTL CKPLVPREHQ FYETLPAEMR RFTPQYKGVV SVRFEEDEDR NLCLIAYPLK GDHGPDIVD NSDCEPKSKL LRWTNKKHHV LETEKSPKDW VRQHRKEEKM KSHKLEEEFE WLKKSEVLYY SVEKKGTVSS QLKHYNPWSM KCHQQQLQRM KENAKHRNQY KFILLNLTC RYEVPCVLDL KMGTRQHGDD ASEEKAANQI RKCQQSTSAV IGVRVCGMQV YQAGTGQLMF MNKYHGRKLS VQGFKEALFQ FFHNGRYLRR ELLGPVLKKL TELKAVLERQ ESYRFYSSSL LVIYDGKEWP EVTLSDAED LEDLSEESAD ESAGAYAYKP LGASSVDVRM IDFAHTTCRL YGEDSVVHEG QDAGYIFGLQ SLIDIVTEIS EESGE
Source	E.coli
Target Names	Ip6k2
Protein Names	Recommended name: Inositol hexakisphosphate kinase 2 Short name= InsP6 kinase 2 EC= 2.7.4.21 Alternative name(s): P(i)-uptake stimulator Short name= PiUS
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein that belongs to the inositol phosphokinase (IPK) family. This protein is likely responsible for the conversion of inositol hexakisphosphate (InsP6) to diphosphoinositol pentakisphosphate (InsP7/PP-InsP5). It may also convert 1,3,4,5,6-pentakisphosphate (InsP5) to PP-InsP4 and affect the growth suppressive and apoptotic activities of interferon-beta in some ovarian cancers. Alternative splicing results in multiple transcript variants encoding different isoforms.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.