



Recombinant Human Tubulin polymerization-promoting protein family member 3 (TPPP3)

Product Code	CSB-EP024117HU
Abbreviation	TPPP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BW30
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AASTD MAGL EESFRKFAIH GDPKASGQEM NGKNWAKLCK DCKVADGKSV TGTDVDIVFS KVKGKSARVI NYEEFKKALE ELATKRFGK SKEEAFDAIC QLVAGKEPAN VGVTKAKTGG AVDRLTDTSR YTGSHKERFD ESGKGKGIAG RQDILDDSGY VSAYKNAGTY DAKVKK
Source	E.coli
Target Names	TPPP3
Protein Names	Recommended name: Tubulin polymerization-promoting protein family member 3 Alternative name(s): TPPP/p20
Expression Region	2-176
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.