



Recombinant Human Golgi phosphoprotein 3-like (GOLPH3L)

Product Code	CSB-BP009669HU
Abbreviation	GOLPH3L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H4A5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTTLTHRARR TEISKNSEKK MESEEDSNWE KSPDNEDSGD SKDIRLTLME EVLLGLKDK EGYTSFWNDC ISSGLRGGIL IELAMRGRIY LEPPTMRKKR LLDRKVLLKS DSPTGDVLLD ETLKHIKATE PTETVQTWIE LLTGETWNPF KLQYQLRNVR ERIAKNLVEK GILTTEKQNF LLFDMTTHPV TNTTEKQRLV KKLQDSVLER WNDPQRMDK RTLALLVLAH SSDVLENVFS SLTDDKYDVA MNRKDLVEL DPEVEGKPS ATEMIWAVLA AFNKS
Source	Baculovirus
Target Names	GOLPH3L
Protein Names	Recommended name: Golgi phosphoprotein 3-like Alternative name(s): GPP34-related protein
Expression Region	1-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. This protein is localized at the Golgi stack and may have a regulatory role in Golgi trafficking.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.