



Recombinant Human Zinc finger protein 395 (ZNF395)

Product Code	CSB-MP026718HU
Abbreviation	ZNF395
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H8N7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASVLSRRLG KRSLGARVL GPSASEGPSA APPSEPLLEG AAPQPFTTSD DTPCQEQPKE VLKAPSTGL QQVAFQPGQK VYVWYGGQEC TGLVEQHSWM EGQVTWVLE QKLQVCCRVE EVWLAELQGP CPQAPPLEPG AQALAYRPVS RNIDVPKRKS DAVEMDEMMA AMVLTSLSCS PVVQSPPGTE ANFSASRAAC DPWKESGDIS DSGSSTTSGH WSGSSGVSTP SPPHPQASPK YLGDAFGSPQ TDHGFETDPD PFLLEPAPR KRKNSVKVMY KCLWPNCVKV LRSIVGIKRH VKALHLGDTV DSDQFKREED FYYTEVQLKE ESAAAAAAAA AGTPVPGTPT SEPAPTPSMT GLPLSALPPP LHKAQSSGPE HPGPESSLPS GALSXSAPGS FWHIQADHAY QALPSFQIPV SPHIYTSVSW AAAPSAACSL SPVRSRSLSF SEPQQPAPAM KSHLIVTSPQ RAQSGARKAR GEAKKCRKVY GIEHRDQWCT ACRWKKACQR FLD
Source	Mammalian cell
Target Names	ZNF395
Protein Names	Recommended name: Zinc finger protein 395 Alternative name(s): HD-regulating factor 2 Short name= HDRF-2 Huntington disease gene regulatory region-binding protein 2 Short name= HD gene regulatory region-binding protein 2 Shor
Expression Region	1-513
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.