



# Recombinant Mouse Suppressor of cytokine signaling 6 (Socs6)

<b>Product Code</b>	CSB-MP022395MO
<b>Abbreviation</b>	Socs6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9JLY0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKKISLKTFR KSFNLSKSKD ETEFMVVQPQ SLAGDFVKDD SLFGSCYGKD MASCDIGSED EKGKNRSKSE SLMGTLKRRL SAKQKTKGKG GTASTDEDTF SSASAPGGLK DVRAPRPIRS TSLRSHHYSPTWPLRPTSS EETCIKMEMR VKALVHAASP GPVNGVRKDL RELQPRELRD LQPEPRPESR CSPSSPGDLS LHLEEHVPVV IGLMSQDYLQ YTVPLDDGMC PLEGPRSCCL DTSSPMEVSA VPLPGASGAF SEDDSHVDQD LVVGPEILVD SSVNNLLIGT TGVMLQSPRG GHDDAPPLSP LLPPMQNNPI QRNFSGLSGP DLHMAESVRC HLNFDPN SAP GVARVYDSVQ SSGPMVVTSL TEELKKLAKQ GWYWGPI TRW EAEGKLANVP DGSFLVRDSS DDRYLLSLSF RSHGKTLHTR IEHSNGRFSF YEQPDVEGHT SIVDLIEHSI RSENGAFCY SRSRLPGSAT YPVRLTNPVS RFMQVRS LQY LCRFVIRQYT RIDLIQKLPL PNKMKDY LQE KHY
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Socs6
<b>Protein Names</b>	Recommended name: Suppressor of cytokine signaling 6 Short name= SOCS-6 Alternative name(s): Cytokine-inducible SH2 protein 4 Short name= CIS-4 Suppressor of cytokine signaling 4 Short name= SOCS-4
<b>Expression Region</b>	1-533
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein contains a SH2 domain and a CIS homolog domain. The protein thus belongs to the cytokine-induced STAT inhibitor (CIS), also known as suppressor of cytokine signaling (SOCS) or STAT-induced STAT inhibitor (SSI), protein family. CIS family members are known to be cytokine-inducible negative regulators of cytokine signaling. The expression of this gene can be induced by GM-CSF and EPO in hematopoietic cells. A high expression level of this gene



was found in factor-independent chronic myelogenous leukemia (CML) and erythroleukemia (HEL) cell lines.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.