



Recombinant Human CREB/ATF bZIP transcription factor (CREBZF)

| | |
|--------------------------|---|
| Product Code | CSB-EP005956HU-B |
| Abbreviation | CREBZF |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9NS37 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MRHSLTKLLA ASGSNSPTRS ESPEPAATCS LPSDLTRAAA GEEETAAAGS PGRKQQFGDE GELEAGRGSR GGVAVRAPSP EEMEEETAIAS LPGEETEDMD FLSGLELADL LDPRQPDWHL DPGLSSPGPL SSSGGGSDSG GLWRGDDDDDE AAAAEMQRFS DLLQRLNGI GGCSSSSSDSG SAEKRRRKSP GGGGGGSGN DNNQAATKSP RKAATAAARL NRLKKKEYVM GLSRVRGLA AENQELRAEN RELGKRVQAL QEESRYLRV LANETGLARL LSRLSGVGLR LTTSLFRDSP AGDHDYALPV GKQKQDLLEE DDSAGGVCLH VDKDKVSVEF CSACARKASS SLKM |
| Source | E.coli |
| Target Names | CREBZF |
| Protein Names | Recommended name: CREB/ATF bZIP transcription factor Alternative name(s): Host cell factor-binding transcription factor Zhangfei Short name= HCF-binding transcription factor Zhangfei |
| Expression Region | 1-354 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.