



Recombinant Human CREB/ATF bZIP transcription factor (CREBZF)

Product Code	CSB-MP005956HU
Abbreviation	CREBZF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NS37
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRHSLTKLLA ASGSNSPTRS ESPEPAATCS LPSDLTRAAA GEEETAAAGS PGRKQQFGDE GELEAGRGSR GGVAVRAPSP EEMEEETAIAS LPGEETEDMD FLSGLELADL LDPRQPDWHL DPGLSSPGPL SSSGGGSDSG GLWRGDDDDDE AAAAEMQRFS DLLQRLNGI GGCSSSSSDSG SAEKRRRKSP GGGGGGSGN DNNQAATKSP RKAATAAARL NRLKKKEYVM GLESVRVGLA AENQELRAEN RELGKRVQAL QEESRYLRV LANETGLARL LSRLSGVGLR LTTSLFRDSP AGDHDYALPV GKQKQDLLEE DDSAGGVCLH VDKDKVSVEF CSACARKASS SLKM
Source	Mammalian cell
Target Names	CREBZF
Protein Names	Recommended name: CREB/ATF bZIP transcription factor Alternative name(s): Host cell factor-binding transcription factor Zhangfei Short name= HCF-binding transcription factor Zhangfei
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.