



Recombinant Human Zinc finger protein 426 (ZNF426)

Product Code	CSB-MP026740HU
Abbreviation	ZNF426
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BUY5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAADLSHGH YLSGDPVCLH EEKTPAGRIV ADCLTDCYQD SVTFDDVAVD FTQEEWTLTD STQRSLYSDV MLENYKNLAT VGGQIIKPSL ISWLEQEESR TVQGGVLQGW EMRLETQWSI LQQDFLRGQT SIGIQLEGKH NGRELCDCEQ CGEVFSEHSC LKTHVRTQST GNTHDCNQYG KDFLTLCEKT STGEKLSEFN QSEKIFSLTP NIVYQRTSTQ EKSFECSHCG KSFINESYLQ AHMRTHNGEK LYEWRNYGPG FIDSTLSVL IETLNAKKPY KCKEKGKGYR YPAYLSIHMR THTGEKPYEC KECGKAFNYS NSFQIHGRTH TGEKPYVCKE CGKAFTQYSG LSMHVRSHSG DKPYECKECG KSFLTSSRLI QHIRTHTGEK PFVCVECGKA FAVSSNLSGH LRTHTEEKAC ECKICGKVFG YPSCLNNHMR THSAQKPYTC KECGKAFNYS THLKIHMRIH TGEKPYECKQ CGKAFSHSSS FQIHERHTHG EKPYECKECG KAFTCSSSFR IHEKTHTEEK PYKCQQCGKA YSHPRSLRRH EQIH
Source	Mammalian cell
Target Names	ZNF426
Protein Names	Recommended name: Zinc finger protein 426
Expression Region	1-554
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.