



Recombinant Chlorocebus aethiops Cathepsin L1 (CTSL1)

Product Code	CSB-BP006193DSU
Abbreviation	CTSL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GKL8
Product Type	Recombinant Protein
Immunogen Species	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
Purity	≥85% (SDS-PAGE)
Sequence	APRSVDWREKGYVTPVKNQGQCGSCWAFSATGALEGQMFRKTGKLVSLSE QNLVDCSGPQGNEGCNGGLMDYAFQYVADNGGLDSEESYPYEATEESCKYN PEYSVANDTGFVDIPKQEKALMKAVATVGPISVAIDAGHESFMFYKEGIYFEPD CSSEDMDHGVLVVGYGFESTESDNSKYWLVKNSWGEEWGMGGYIKMAKDR RNHCGIASAASYPTV
Source	Baculovirus
Target Names	CTSL
Protein Names	Recommended name: Cathepsin L1 EC= 3.4.22.15 Alternative name(s): Major excreted protein Short name= MEPCleaved into the following 2 chains: 1. Cathepsin L1 heavy chain 2. Cathepsin L1 light chain
Expression Region	114-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.