



Recombinant Human Zinc finger protein 350 (ZNF350)

Product Code	CSB-BP026694HU
Abbreviation	ZNF350
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZX5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MIQAQESITLEDVAVDFTWEEWQLLGAAQKDLYRDVMLENYSNLVAVGYQAS KPDALFKL EQGEQLWTIEDGIHSGACSDIWKVDHVLRLQSESLVNRKPCHEHDAFENIV HCSKSQF LLGQNHDIFDLRGKSLKSNLTLVNQSKGYEIKNSVEFTGNGDSFLHANHERLH TAIKFP SQKLSTKSQFISPKHQKTRKLEKHHVCSECGKAFIKKSWLTDHQMHTGKPK HRCSLCE KAFSRKFMLTEHQRTHTGKPYECPECGKAFLLKSRNLNIHQKTHHTGKPYICS ECGKGF QKGNLIVHQRIHTGKPYICNECGKGFQKTCIAHQRFHTGKTPFVCSECGKS CSQKSG LIKHQRIHTGKPFECSECGKAFSTKQKLIVHQRTHTGERPYGCNECGKAFAY MSCLVKH KRIHTREKQEAQVENPPAERHSSLHTSDVMQEKNSANGATTQVPSVAPQTS LNISGLLA NRNVVLVGQPVVRCAASGDNRGFAQDRNLVNAVNVVPSVINYVLFYVTENP</p>
Source	Baculovirus
Target Names	ZNF350
Protein Names	Recommended name: Zinc finger protein 350 Alternative name(s): KRAB zinc finger protein ZFQR Zinc finger and BRCA1-interacting protein with a KRAB domain 1 Zinc finger protein ZBRK1
Expression Region	1-532
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.