



# Recombinant Human Zinc finger protein 350 (ZNF350)

<b>Product Code</b>	CSB-EP026694HU-B
<b>Abbreviation</b>	ZNF350
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9GZX5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MIQAQESITLEDVAVDFTWEEWQLLGAAQKDLYRDVMLENYSNLVAVGYQAS            KPDALFKL            EQGEQLWTIEDGIHSGACSDIWKVDHVLRLQSESLVNRKPCHEHDAFENIV            HCSKSQF            LLGQNHDFDLRGKSLKSNLTLVNQSKGYEIKNSVEFTGNGDSFLHANHERLH            TAIKFP            SQKLSTKSQFISPKHQKTRKLEKHHVCSECGKAFIKKSWLTDHQMHTGEKP            HRCSLCE            KAFSRKFMLTEHQRTHTGEKPYECPECGKAFLLKSRNLNIHQKTHHTGEKPYICS            ECGKGF            QKGNLIVHQRIHTGEKPYICNECGKGFQKTCIAHQRFHTGKTPFVCSECGKS            CSQKSG            LIKHQRIHTGEKPFECSECGKAFSTKQKLIVHQRTHTGERPYGCNECGKAFAY            MSCLVKH            KRIHTREKQEAACKVENPPAERHSSLHTSDVMQEKNSANGATTQVPSVAPQTS            LNISGLLA            NRNVVLVGQPVVRCAASGDNRGFAQDRNLVNAVNVVPSVINYVLFYVTENP</p>
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF350
<b>Protein Names</b>	Recommended name: Zinc finger protein 350 Alternative name(s): KRAB zinc finger protein ZFQR Zinc finger and BRCA1-interacting protein with a KRAB domain 1 Zinc finger protein ZBRK1
<b>Expression Region</b>	1-532
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.