



# Recombinant Human Zinc finger protein 703 (ZNF703)

<b>Product Code</b>	CSB-MP026966HU
<b>Abbreviation</b>	ZNF703
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H7S9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SDSPAGSNP RTPESSGSGS GGGGKRPAVP AAVSLLPPAD PLRQANRLPI RVLKMLSAHT GHLLHPEYLQ PLSSTPVSPI ELDAKKSPLA LLAQTCSQIG KPDPSSSKL NSVAAAANGL GAEKDPGRSA PGAASAAAAL KQLGDSPAED KSSFKPYSKG SGGGDSRKDS GSSSVSSTSS SSSSSPGDKA GFRVPSAACP PFPPHGAPVS ASSSSSSPGG SRGGSPHHSD CKNGGGVGGG ELDKKDQEPK PSPEPAAVSR GGGGEPGAHG GAESGASGRK SEPPSALVGA GHVAPVSPYK PGHSVFPLPP SSIGYHGSIV GAYAGYPSQF VPGLDPSKSG LVGGQLSGGL GLPPGKPPSS SPLTGASPPS FLQGLCRDPY CLGGYHGASH LGGSSCSTCS AHDPAGPSLK AGGYPLVYPG HPLQPAALSS SAAQAALPGH PLYTYGFMLQ NEPLPHSCNW VAASGPCDKR FATSEELLSH LRTHTALPGA EKLLAAYPGA SGLGSAAAAA AAAASCHLHL PPPAAPGSPG SLSLRNPHTL GLSRYHPYGK SHLSTAGGLA VPSLPTAGPY YSPYALYGQR LASASALGYQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF703
<b>Protein Names</b>	Recommended name: Zinc finger protein 703 Alternative name(s): Zinc finger elbow-related proline domain protein 1
<b>Expression Region</b>	2-590
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.