



Recombinant Human Zinc finger protein 703 (ZNF703)

Product Code	CSB-MP026966HU
Abbreviation	ZNF703
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H7S9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SDSPAGSNP RTPESSGSGS GGGGKRPAVP AAVSLLPPAD PLRQANRLPI RVLKMLSAHT GHLLHPEYLQ PLSSTPVSPI ELDAKKSPLA LLAQTCSQIG KPDPPSSKL NSVAAAANGL GAEKDPGRSA PGAASAAAAL KQLGDSPAED KSSFKPYSKG SGGGDSRKDS GSSSVSSTSS SSSSSPGDKA GFRVPSAACP PFPPHGAPVS ASSSSSSPGG SRGGSPHSD CKNGGGVGGG ELDKKDQEPK PSPEPAAVSR GGGGEPGAHG GAESGASGRK SEPPSALVGA GHVAPVSPYK PGHSVFLPP SSIGYHGSIV GAYAGYPSQF VPGLDPSKSG LVGGQLSGGL GLPPGKPPSS SPLTGASPPS FLQGLCRDPY CLGGYHGASH LGGSSCSTCS AHDPAGPSLK AGGYPLVYPG HPLQPAALSS SAAQAALPGH PLYTYGFMLQ NEPLPHSCNW VAASGPCDKR FATSEELLSH LRHTALPGA EKLLAAYPGA SGLGSAAAAA AAAASCHLHL PPPAAPGSPG SLSLRNPHL GLSRYHPYGK SHLSTAGGLA VPSLPTAGPY YSPYALYGQR LASASALGYQ
Source	Mammalian cell
Target Names	ZNF703
Protein Names	Recommended name: Zinc finger protein 703 Alternative name(s): Zinc finger elbow-related proline domain protein 1
Expression Region	2-590
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.