



Recombinant Human Integrator complex subunit 10 (INTS10)

Product Code	CSB-BP011758HU
Abbreviation	INTS10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NVR2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSAQQGDCEFL VQRARELVPQ DLWAAKAWLI TARSLYPADF NIQYEMYTIE RNAERTATAG RLLYDMFVNF PDQPVVWREI SIITSALRND SQDKQTQFLR SLFETLPGRV QCCEMLLKVTE QCFNTLERSE MLLLLLRRFP ETVVQHGVGL GEALLEAETI EEQESPVNCF RKLFCVDVLP LIINNHVRL PANLLYKYLN KAAEFYINYV TRSTQIENQH QGAQDTS DLM SPSKRSSQKY IIEGLTEKSS QIVDPWERLF KILNVVGMRC EWQMDKGRRS YGDILHRMKD LCRYMNNFDS EAHAKYKNQV VYSTMLVFFK NAFQYVNSIQ PSLFQGP NAP SQVPLVLLED VSNVYGDVEI DRNKHIHKKR KLAEGREKTM SSDDEDCSAK GRNRHIVV NK AELANSTEV L ESFKLARESW ELLYSLEFLD KEFTRICLAW KTDTWLWLRI FLTDMIYQG QYKKAIASLH HLAALQGSIS QPQITGQGTL EHQRALIQLA TCHFALGEYR MTCEKVL DLM CYMVLPIQDG GKSQEEPSKV KPKFRKGSDL KLLPCTSKAI MPYCLHLMLA CFKLR AFTDN RDDMALGHVI VLLQQEWPRG ENLFLKAVNK ICQQGNFQYE NFFNYVTNID MLEEFAYLRT QEGGKIHLEL LPNQGMLIKH HTVTRGITKG VKEDFRLAME RQVSRGENL MVVLHRFCIN EKILLQTLT
Source	Baculovirus
Target Names	INTS10
Protein Names	Recommended name: Integrator complex subunit 10 Short name= Int10
Expression Region	1-710
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.