



# Recombinant Mouse Transcription factor EB (Tfeb)

<b>Product Code</b>	CSB-BP023429MO
<b>Abbreviation</b>	Tfeb
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9R210
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MASRIGLRMQ LMREQAQQEE QRERMQQQAV MHYMQQQQQQ QQQLGGPPTP AINTPVHFQS PPPVPGEVLK VQSYLENPTS YHLQQSQHQK VREYLSETYG NKFAAHVSPA QGSPKPAPAA SPGVRAGHVL STSAGNSAPN SPMAMLISS NPEKEFDDVI DNIMRLDSVL GYINPEMQMP NTLPLSSSHL NVYSGDPQVT ASMVGVTSSS CPADLTQKRE LTDAESRALA KERQKKNHN LIERRRRFNI NDRIKELGML IPKANDLDVR WNKGTILKAS VDYIRRMQKD LQKSRELENH SRRLEMTNKQ LWLRIQELM QARVHGLPTT SPSGVNMAEL AQQVVKQELP SEDGPGREALM LGPEVPEPEQ MPALPPQAPL PSAAQPQSPF HHLDFSHGLS FGGGGDEGPT GYPDTLGTTEH GSPFPNLSKK DLDLMLLDDS LLPLASDPLF STMSPEASKA SSRSSFSME EGDVL
<b>Source</b>	Baculovirus
<b>Target Names</b>	Tfeb
<b>Protein Names</b>	Recommended name: Transcription factor EB
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.