



Recombinant Human Nucleolar protein 10 (NOL10)

Product Code	CSB-EP015918HU-B
Abbreviation	NOL10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSC4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MQVSSLNEVK IYSLSCGKSL PEWLSDRKKR ALQKKDVDVR RRIELIQDFE MPTVCTTIKV SKDGQYILAT GTYKPRVRCY DTYQLSLKFE RCLDSEVVTF EILSDDYSKI VFLHNDRYIE FHSQSGFYK TRIPKFRGDF SYHYPSCDLY FVGASSEVYR LNLEQGRYLN PLQTDAEENN VCDINSVHGL FATGTIEGRV ECWDPRTNR VGLLDCALNS VTADSEINSL PTISALKFNG ALTMAVGTTT GQVLLYDLRS DKPLLVKDHQ YGLPIKSVHF QDSLILLSA DSRIVKMWNK NSGKIFTSLE PEHDLNDVCL YPNSGMLLTA NETPKMGIYY IPVLGPAPRW CSFLDNLTEE LEENPESTVY DDYKFVTKKD LENLGLTHLI GSPFLRAYMH GFFMDIRLYH KVKLMVNPFA YEEYRKDKIR QKIEETRAQR VQLKKLPKVN KELALKLIEE EEEKQKSTWK KVKSLPNIL TDDRFKVMFE NPFDQVDEES EEFRLNPLV SKISEKRKKK LRLLEQQELR EKEEEEEPEG KPSDAESSES SDDEKAWVEE VRKQRLLQQ EEKVKRQERL KEDQQTVLKP QFYEIKAGEE FRSFKDSATK QKLMNKTLED RLKIEAKNGT LSVSDTTVGS KQLTFTLKRS EQQKQQAEE KLHRQERKRL RRSAGHLKSR HKRGRSFH
Source	E.coli
Target Names	NOL10
Protein Names	Recommended name: Nucleolar protein 10
Expression Region	1-688
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.