



Recombinant Human Zinc finger protein 331 (ZNF331)

Product Code	CSB-MP026680HU
Abbreviation	ZNF331
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NQX6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAQGLVTFAD VAIDFSQEEW ACLNSAQRDL YWDVMLENYS NLVSLDLESA YENKSLPTEK NIHEIRASKR NSDRRSKSLG RNWICEGTLE RPQRSRGRYV NQMIINYVKR PATREGTPPR THQRHHKENS FECKDCGKAF SRGYQLSQHQ KIHTGEKPYE CKECKKAFRW GNQLTQHQQI HTGEKPYECK DCGKAFRWGS SLVIHKRIHT GEKPYECKDC GKAFRRGDEL TQHQRFHTGE KDYECKDCGK TFSRVYKLIQ HKRIHSGEKP YECKDCGKAF ICGSSLIQHK RIHTGEKPYE CQECGKAFTR VNYLTQHQQI HTGEKPHECK ECGKAFRWGS SLVKHERIHT GEKPYKCTEC GKAFNCGYHL TQHERIHTGE TPYKCKEKGK AFIYGSSLVK HERIHTGVKP YGCTECGKSF SHGHQLTQHQQI KTHSGAKSYE CKECGKACNH LNHLREHQRI HNS
Source	Mammalian cell
Target Names	ZNF331
Protein Names	Recommended name: Zinc finger protein 331 Alternative name(s): C2H2-like zinc finger protein rearranged in thyroid adenomas Zinc finger protein 361 Zinc finger protein 463
Expression Region	1-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.