



Recombinant Human REST corepressor 3 (RCOR3)

Product Code	CSB-EP019517HU-B
Abbreviation	RCOR3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P2K3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRVGAEYQAR IPEFDPGATK YTDKDNGGML VWSPYHSIPD AKLDEYIAIA KEKHGYNVEQ ALGMLFWHKH NIEKSLADLP NFTPFPDEWT VEDKVLFEQA FSFHGKSFHR IQQMLPDKTI ASLVKYYYSW KKTRSRTSLM DRQARKLANR HNQGDSDDDDV EETHPMDGND SDYDPKKEAK KEGNTEQPVQ TSKIGLGRRE YQSLQHRHHS QRSKCRPPKG MYLTQEDVVA VSCSPNAANT ILRQLDMELI SLKRQVQNAK QVNSALKQKM EGGIEEFKPP ESNQKINARW TTEEQLLAVQ GVRKYGKDFQ AIADVIGNKT VGQVKNFFVN YRRRFNLEEV LQEWAEQGT QASNGDASTL GEETKSASNV PSGKSTDEEE EAQTPQAPRT LGPSPAPSS TPTPTAPIAT LNQPPLLRP TLPAAPALHR QPPPLQQQAR FIQPRPTLNQ PPPPLIRPAN SMPRLNRP VLSTVGGQQP PSLIGIQTDS QSSLH
Source	E.coli
Target Names	RCOR3
Protein Names	Recommended name: REST corepressor 3
Expression Region	1-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.