



Recombinant Mouse Prolyl endopeptidase (Prep)

Product Code	CSB-MP018664MO
Abbreviation	Prep
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9QUR6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MLSFQYPDVY RDETSVQEYH GHKICDPYSW LEDPDSEQTK AFVEAQNKIT VPFLEQCPIR GLYKERMTEL YDYPKYSCHF KKGKRYFYFY NTGLQNQRVL YVQDSLEGEA RVFLDPNTLS DDGTVALRGY AFSEDGEYFA YGLSASGSDW VTIKFMKVDG AKELPDVLER VKFTCMAWTH DGKGMFYNSY PQQDGKSDGT ETSTNLHQKL CYHVLGTDQS EDILCAEFPD EPKWMGGAEL SDDGRYVLLS IWEGCDPVNR LWYCDLQQEP NGITGILKWV KLIDNFEGEY DYVTNEGTVF TFKTNRNSPN YRLINIDFTD PDESKWKVLV PEHEKDVLEW VACVRSNFLV LCYLHDVKNI LQLHDLTGGA LLKTFPLDVG SVVGYSGRKK DSEIFYQFTS FLSPGVYHC DLTKEELEPM VFREVTVKGI DAADYQTIQI FYPSKDGTKI PMFIVHKKGI KLDGSHPAFL YGYGGFNISI TPNYSVSRLI FVRHMGGVLA VANIRGGGEY GETWHKGGIL ANKQNCFFDDF QCAAAYLIKE GYTSPKRLTI NGGSNGGLLV AACANQRPDL FGCVIAQVGV MDMLKFHKFT IGHAWTTDYG CSDTKQHFEEW LLKYSPLHNV KLPEADDIQY PSMLLLTADH DDRVVPLHSL KFIATLQYIV GRSRKQSNPL LIHVDTKAGH GAGKPTAKVI EEVSDMFAFI ARCLNIEWIQ
Source	Mammalian cell
Target Names	Prep
Protein Names	Recommended name: Prolyl endopeptidase Short name= PE EC= 3.4.21.26 Alternative name(s): Post-proline cleaving enzyme
Expression Region	1-710
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a cytosolic prolyl endopeptidase that cleaves peptide bonds on the C-terminal side of prolyl residues within peptides that are up to approximately 30 amino acids long. Prolyl endopeptidases have been reported to be involved in the maturation and degradation of peptide hormones and neuropeptides.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.