



Recombinant Human Zinc finger protein 282 (ZNF282)

Product Code	CSB-YP026638HU
Abbreviation	ZNF282
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UDV7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MQFVSTRPQP QQLGIQGLGL DSGSWSWAQA LPPEEVCHQE PALRGEMAEG MPPMQAQEWD MDARRPMPFQ FPPFPDRAPV FPDRMMREPQ LPTAEISLWT VVAAIQAVER KVDAQASQLL NLEGRTGTAE KKLADCEKTA VEFGNHMESK WAVLGTLLQE YGLLQRRLLEN LENLLRNRNF WVLRLPPGSK GEAPKVPVTF VDIAYVFSER EWKNLDEWQK ELYNNLVKEN YKTLMSLDAE GSVPKPDAPV QAEPREEPCV WEQRHPEERE IPMDPEAGAE PLVPAQDASS QVKREDTLCV RGQRGLEERA IPTESITDSP ISAQDLLSRI KQEEHQCVWD QQDLADRIP TDPNSESIS AHDILSWIKQ EEQPYPWGPR DSMDGELGLD SGPSDSLMLV KNPPPAPPQP QPQPQPPQPQ LQSQPQPQSL PPIAVAENPG GPPSRGLLDD GFQVLPGERG SGEAPPGGDR STGGGGGGDGG GGGGGAEAGT GAGGGCGSCC PGGLRRSLLL HGARSKPYSC PECGKSFGVR KSLIIHHRSH TKERPYECAE CEKSFNCHSG LIRHQMTHRG ERPYKCSECE KTYSRKEHLQ NHQRLHTGER PFQCALCGKS FIRKQNLLKH QRIHTGERPY TCGECGKSFR YKESLKDHLR VHSGGPGPGA PRQLPPPPER D </p>
Source	Yeast
Target Names	ZNF282
Protein Names	Recommended name: Zinc finger protein 282 Alternative name(s): HTLV-I U5RE-binding protein 1 Short name= HUB-1
Expression Region	1-671
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.