



Recombinant Human INO80 complex subunit B (INO80B)

Product Code	CSB-MP011725HU
Abbreviation	INO80B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C086
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSKLWRRGST SGAMEAPEPG EALELSLAGA HGHGVHKKKH KKHKKKHKKK HHQEEDAGPT QPSPAKPQLK LKIKLGGQVL GTKSVPTFTV IPEGPRSPSP LMVVDNEEEP MEGVPLEQYR AWLDEDSNLS PSPLRDLSGG LGGQEEEEEQ RWLDALEKGE LDDNGDLKKE INERLLTARQ RALLQKARSQ PSPMLPLPVA EGCPPPALTE EMLLKREERA RKRRQLQAARR AEEHKNQTIE RLTKTAATSG RGGRGGARGE RRGGRAAAPA PMVRYCSGAQ GSTLSFPPGV PAPTAVSQRP SPSGPPPRCS VPGCPHPRRY ACSRTGQALC SLQCYRINLQ MRLGGPEGPG SPLLAT
Source	Mammalian cell
Target Names	INO80B
Protein Names	Recommended name: INO80 complex subunit B Alternative name(s): High mobility group AT-hook 1-like 4 IES2 homolog Short name= hles2 PAP-1-associated protein 1 Short name= PAPA-1 Zinc finger HIT domain-containing protein 4
Expression Region	1-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.