



# Recombinant Rat Mesothelin (Msln)

<b>Product Code</b>	CSB-EP015044RA-B
<b>Abbreviation</b>	Msln
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9ERA7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DTEQKACPPGKEPNVVDENLIFYQNWELEACVDGTLLAGQMDLVNEIPFTYEQ LSIFKHKLDKTYPQGYPESLIKQLGHFFRYVSPEDIRQWNVTSPDTVNTLLKVS KGQKMDAQVIALVACYLRGGGKLEDEDIVKALDNIPLSYLCDFSPQDLHAIPSSV MWLVGLHDLDKCSQRHLGILYQKACSAFQNVSGLEYFEKIRTFLLGGASREDLR ALSQHNVSMDIATFKKLQVDALVGLSVAEVQKLLGPHIGDLKTEEDKSPVRDW LFRQQQKDLDSLGLGLQGIPNGYLILDFNVREAFS
<b>Source</b>	E.coli
<b>Target Names</b>	Msln
<b>Protein Names</b>	Recommended name: Mesothelin Alternative name(s): Pre-pro-megakaryocyte-potentiating factor Protein expressed in renal carcinoma Cleaved into the following 2 chains: 1. Megakaryocyte-potentiating factor Short name= 2. MPF 3. Mes
<b>Expression Region</b>	298-600
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	An antibody that reacts with ovarian cancers and mesotheliomas was used to isolate a cell surface antigen named mesothelin. Although the function of mesothelin is unknown, it may play a role in cellular adhesion and is present on mesothelium, mesotheliomas, and ovarian cancers. This gene encodes a preproprotein that is cleaved into two products, megakaryocyte potentiating factor and mesothelin. Alternative splicing results in multiple transcript variants.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.