



Recombinant Human Docking protein 5 (DOK5)

Product Code	CSB-BP007111HU
Abbreviation	DOK5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9P104
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASNFNDIVK QGYVRIRSRR LGIYQRCWL V FKKASSKGP K RLEKFSDERA AYFRCYHKVT ELNNVKNVAR LPKSTKKHAI GIYFNDDTSK TFACESDLEA DEWCKVLQME CVGTRINDIS LGEPDLLATG VEREQSERFN VYLMPSPNLD VHGEALQIT YEYICLWDVQ NPRVKLISWP LSALRRYGRD TTWFTFEAGR MCETGGLFI FQTRDGEAIY QKVHSAALAI AEQHERLLQS VKNSMLQMKM SERAASLSTM VPLPRSAYWQ HITRQHSTGQ LYRLQDVSSP LKLHRTETFP AYRSEH
Source	Baculovirus
Target Names	DOK5
Protein Names	Recommended name: Docking protein 5 Alternative name(s): Downstream of tyrosine kinase 5 Insulin receptor substrate 6 Short name= IRS-6 Short name= IRS6
Expression Region	1-306
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the DOK family of membrane proteins, which are adapter proteins involved in signal transduction. The encoded protein interacts with phosphorylated receptor tyrosine kinases to mediate neurite outgrowth and activation of the MAP kinase pathway. In contrast to other DOK family proteins, this protein does not interact with RASGAP.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.