



Recombinant Human Leucine carboxyl methyltransferase 1 (LCMT1)

Product Code	CSB-BP012808HU
Abbreviation	LCMT1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UIC8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MATRQRESSI TSCCSTSSCD ADDEGVRGTC EDASLCKRFA VSIGYWHDPY IQHFVRLSKE RKAPEINRGY FARVHGVSQI IKAFLRKTEC HCQIVNLGAG MDTTFWRLKD EDLLPSKYFE VDFPMIVTRK LHSIKCKPPL SSPILELHSE DTLQMDGHIL DSKRYAVIGA DLRDLSELEE KLKKNMNTQ LPTLLIAECV LVYMTPEQSA NLLKWAANSF ERAMFINYEQ VNMGDRFGQI MIENLRRRQC DLAGVETCKS LESQKERLLS NGWETASAVD MMELYNRLPR AEVSRIESLE FLDEMELLEQ LMRHYCLCWA TKGGNELGLK EITY
Source	Baculovirus
Target Names	LCMT1
Protein Names	Recommended name: Leucine carboxyl methyltransferase 1 EC= 2.1.1.233 Alternative name(s): Protein-leucine O-methyltransferase [Phosphatase 2A protein]-leucine-carboxy methyltransferase 1
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.