



Recombinant Human Zinc finger protein 230 (ZNF230)

Product Code	CSB-BP026598HU
Abbreviation	ZNF230
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UIE0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTTFKEAVTF KDVAVFFTEE ELGLLDPAQR KLYQDVMLEN FTNLLSVGHQ PFHPPHFLRE EKFWMMETAT QREGNSGGKT IAEAGPHEDC PCQQIWEQTA SDLTQSQDSI INNSHFFEQQ DVPSQVEAGL SIIHTGQKPS QNGKCKQSFS DVAIFDPPQQ FHSGEKSHTC NECGKSFCYI SALRIHQRVH LREKLSKCDM RGKEFSQSSC LQTRERVHTG EKPFKCEQCG KGFRRCRAILQ VHCKLHTGEK PYICEKCGRA FIHDFQLQKH QIIHTGEKPF KCEICGKSFC LRSSLNRHCM VHTAEKLYKS EECGKGFTDS LDLHKHQIIH TGQKPYNCKE CGKSFRWSSY LLIHQRIHSG EKPYPCEECEG KGYISKGLN LHQRVHTGER PYNCKEKGKS FSRASSILNH KKLHCRKKPF KCEDCGKRLV HRSFCKDQQG DHNGENSSKC EDCGKRYKRR LNLDIILSLF LNDM
Source	Baculovirus
Target Names	ZNF230
Protein Names	Recommended name: Zinc finger protein 230 Alternative name(s): Zinc finger protein FDZF2
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.